


FOREHEAD THERMOMETER



- ✓ NON-CONTACT MEASUREMENT, JUST PRESS AND READ
- ✓ QUICK RESPONSE, LESS THAN 1 SECOND
- ✓ SANITARY, SAFE AND ACCURATE
- ✓ AUTO POWER OFF SAVES BATTERY
- ✓ AVAILABLE IN °F OR °C



FEATURES:

- Non contact thermometer for monitoring forehead temperature
- A great portable device to prevent the spread of **SARS** and other viruses
- Safe and sanitary
- Auto power off after 4 seconds of inactivity to preserve battery
- Easy to operate
- Fast response
- High accuracy
- Preset to °F or °C

SPECIFICATIONS:

- **IRT77F** temperature range: 86° to 122°F
- **IRT77C** temperature range: 30° to 50°C
- Accuracy: ±1°F (±0.5°C)
- Response time: 500 m/sec
- Resolution: 0.1°
- Dimensions: 6.7" x 1.7" x 1.6"
- Power source: 1 "9-volt" battery



Carrying Case (optional)

* Note: This thermochecker is only for monitoring forehead temperature and should not be used to judge a fever.

LEADERS IN TEMPERATURE, HUMIDITY & AIR FLOW TECHNOLOGY