

Specifications TGA L87 MSB

	MSB-TGA-2400	MSB-TGA-1600	MSB-TGA-1800	MSB-TGA-1000	MSB-TGA-1500
Resolution	10 µg / 1 µg	10 µg / 1 µg	10 µg / 1 µg	10 µg / 1 µg	10 µg / 1 µg
Range	25 g / 10 g	25 g / 10 g	25 g / 10 g	25 g / 10 g	25 g / 10 g
Vacuum		10 ⁻⁵ mbar*			
Pressure max.	1 bar	5 bar	50 bar	150 bar	150 bar
Temperature range	-150 °C up to 2400 °C**	-150 °C up to 1000 °C Optional: RT up to 1600 °C**	RT up to 1800 °C	RT up to 1000 °C	RT up to 1500 °C
Dosing system	static & dynamic	static & dynamic dynamic	static & dynamic	static & dynamic	static & dynamic
Options	Vapor	Vapor	Vapor	Vapor	Vapor

* Level can be controlled

** Depends on furnace