



pushing boundaries

Linseis Messgeräte GmbH | Vielitzerstr. 43 | D-95100 Selb

Specifications DSC L63

	DSC L92
Temperature range	-40°C to 160°C
Scanning rate	0.001 to 5°C/min
Temperature precision regulation	100 µ°C
Measurement time	30 to 60 min
Sensor type	Peltier Element
Unique Features	Direct T sample measurement Joule effect calibration
Measurement range	±150 mW
RMS noise	0.05 µWa
Equipment dimensions	L*w*h = 900*700*500 mm