



pushing boundaries

Linseis Messgeräte GmbH | Vielitzerstr. 43 | D-95100 Selb

	DIL L73
Temperature range	-180 °C up to 500 °C/ -180 °C up to 1000 °C/ RT up to 1000 °C
Resolution	0.3 nm*
Heating/Cooling rates	0.01 K/min to 50 K/min**
Sample length	up to 50 mm
Sample diameter	up to 7 mm
Sample holder	fused silica
Atmosphere	inert, oxic., red., vac.
Interfaces	USB

* tested in laboratory environments

** depends on furnace

UST-Id.Nr.DE 132 951 806 | TELEFON +49 9287 880-0 | FAX +49 9287 70488 | E-MAIL: info@linseis.de |

Bankkonto: Sparkasse Hochfranken IBAN DE 89 7805 0000 0760 0044 65 BIC BYLADEM1HOF

Bankkonto: Commerzbank IBAN DE 39 7004 0041 0228 8256 00 BIC COBADEFFXXX

Geschäftsleitung: Dipl. Phys. Claus Linseis • Registergericht: HRB HOF ABT. B NR: 299