



pushing boundaries

Linseis Messgeräte GmbH | Vielitzerstr. 43 | D-95100 Selb

Specifications UDSC L64-LT

	UDSC L64-LT
Temperature range	From -170°C to 50 °C
Ramp speed	0,001 to 3 °C/min
Control accuracy	100µ°C
Sample volume	10 to 100 µL
Sensitivity	At 0 °C, 700 µV/mW At -150 °C, 300 µV/mW
Unique features	High sensitive patented sensor
Equipment dimensions	L*w*h = 710*650*550 mm