



pushing boundaries

Linseis Messgeräte GmbH | Vielitzerstr. 43 | D-95100 Selb

Specifications CAL L92

	CAL L92
Temperature range	-60°C* to 170°C
Scanning rate	0.001 to 5°C/min
Temperature precision regulation	100 µ°C
Measurement time	30 to 60 min
Sensor type	Peltier Element
Sample volume	up to 850 mL **
Unique Features	Direct T sample measurement Joule effect calibration
Measurement range	±150 mW
RMS noise	0.05 µW
Equipment dimensions	L*w*h = 700*650*500 mm

* Depends on chiller

** Depends on vessel type